

A simple and integrated workflow for deep proteomic and transcriptomic analysis of sorted cell populations

Authors: Kit Fuhrman¹, Gary DosSantos², Gokhan Demirkan¹, Gary Geiss¹, Doug Hinerfeld¹

¹NanoString Technologies, Seattle, WA

²Reachbio Research Labs, Spokane, WA

In addition to the long standing efforts to understand and manipulate the immune system in the treatment of autoimmune and other immune-related diseases, it is increasingly becoming a direct target for cancer therapeutics. The immense heterogeneity in immune cell phenotype and function and corresponding role in health and disease requires increasingly sophisticated analytical approaches to analysis. Recent advances in flow and mass cytometry have greatly expanded the number of immune cell parameters that can be interrogated, resulting in an improved understanding of the immune system heterogeneity. These technologies, however, remain limited in the number and types of analytes that can be examined in a single clinical sample. Expanding these parameters will enable the discovery of novel biology, potentially leading to the discovery of new therapeutic targets and biomarker signatures.

The NanoString nCounter® platform enables the highly multiplexed digital detection of both RNA and protein from a single biological specimen. Specifically, the nCounter® Vantage 3D™ RNA:Protein Immune Cell Profiling Assay interrogates 30 cell surface proteins and 770 immune-related RNA starting with as few as 20,000 cells in suspension. We demonstrate here the development of a streamlined workflow that seamlessly integrates standard immune cell flow sorting with downstream nCounter® analysis (Figure 1).

By co-staining PBMC with both fluorescently-labeled and DNA-barcoded antibodies, CD8+, CD3+, CD4+ T cells, and CD19+ B cells were isolated by flow cytometry followed by analysis of additional proteins and 770 RNA from each sorted population on the nCounter platform. Demonstrating the value of this workflow in analyzing rare cell populations, the number of target cells were titrated to determine the sensitivity of the workflow (Figure 2). Without the requirement for additional molecular biology methods, such as RNA purification or sequencing library construction, this method is ideally suited for incorporation into any cell sorting workflow.

Sorted Cell Analysis




FIGURE 1 Workflow for the isolation and deep proteomic and genomic analysis of defined cell sorted populations.

High Parameter Protein Analysis from Low Cell Inputs




FIGURE 2 Correlation of the expression of 30 protein markers in defined immune cell populations across different low cell inputs. Values on X and Y axis are normalized counts.

References

1. Geiss, GK, Bumgarner, RE, et al. (2008) Direct multiplexed measurement of gene expression with color-coded probe pairs. *Nat Biotechnol.* Mar;26(3):317-25.
2. Ullal, AV, Peterson, V, et al. (2014) Cancer cell profiling by barcoding allows multiplexed protein analysis in fine-needle aspirates. *Sci Transl Med.* 2014 Jan 15;6(219)
3. Abey, SK, Yuana, Y, et al. (2016) Lysozyme association with circulating RNA, extracellular vesicles, and chronic stress. *BBA Clin.* 2016 Dec 20;7:23-35.

To learn more about 3D Biology™ Technology visit 3d.nanostring.com

NanoString Technologies, Inc.

530 Fairview Avenue North
Seattle, Washington 98109

CONTACT US

info@nanostring.com
Tel: (888) 358-6266
Fax: (206) 378-6288
www.nanostring.com

SALES CONTACTS

United States: us.sales@nanostring.com
EMEA: europe.sales@nanostring.com
Asia Pacific & Japan: apac.sales@nanostring.com
Other Regions: info@nanostring.com

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

©2017 NanoString Technologies, Inc. All rights reserved. NanoString, the NanoString logo, Vantage 3D, 3D Biology and nCounter are trademarks of NanoString Technologies, Inc. in the US and certain other countries.